
SPM MEDICARE PVT. LTD.
(AN EN ISO 13485 Certified Company)
Noida, Uttar Pradesh-201305, India
Customer Care No.: +91-8407070718
E-mail: info@spmmedicare.com
Website: www.spmmedicare.com

MADE IN INDIA


SPMDHF/23/PD/125/Annex-1UR/EU/Rev 01

EASY-FLUSH[®] **Flush Syringe**

Prefilled 0.9% Normal Saline Solution

EN	INSTRUCTIONS FOR USE
DE	GEBRAUCHSANWEISUNG
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
BS	UPUTSTVO ZA UPOTREBU
CS	NÁVOD K POUŽITÍ
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES	INSTRUCCIONES DE USO
ET	KASUTUSJUHEND
FI	KÄYTTÖOHJEET
FR	MODE D'EMPLOI
HR	UPUTE ZA UPOTREBU
HU	HASZNÁLATI UTASÍTÁS
IT	ISTRUZIONI PER L'USO
LT	NAUDOJIMO INSTRUKCIJOS
LV	LIETOŠANAS INSTRUKCIJA
NL	GEbruiksaanwijzing
NO	BRUKSANVISNING
PL	INSTRUKCJA OBSŁUGI
PT	INSTRUÇÕES DE USO
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
SL	NAVODILA ZA UPORABO
SR	UPUTSTVO ZA UPOTREBU
SV	BRUKSANVISNING
TR	KULLANIM TALIMATLARI
UK	ІНСТРУКЦІЯ З ВИКОРИСТАННЯ

EN INSTRUCTIONS FOR USE

 **EASY-FLUSH® FLUSH SYRINGE** (PRE-FILLED, CONTAINS 0.9% SODIUM CHLORIDE SOLUTION)

DESCRIPTION AND USE:

The pre-filled Easy-Flush flushing syringe is a ready-to-use medical device available in various volumes. It is made of high-quality polypropylene and contains a sterile, pyrogen-free 0.9% isotonic sodium chloride solution (NaCl). The syringe's contents are sterile, non-toxic, and pyrogen-free. It is intended solely for flushing compatible intravenous lines and access devices. It is not designed for reconstituting dry products or diluting medications as part of intravenous sodium chloride therapy.

CONTRAINDICATIONS:

Do not use if the patient has a known allergy to any of the components, materials, or the 0.9% sodium chloride solution, as this may cause allergic reactions, including anaphylaxis. Do not use if intravenous administration systems or catheters are contaminated with microorganisms or if the skin around the catheter site is inflamed. Do not use for patients with hypernatremia or fluid retention, as the administration of sodium or chloride may be harmful to health.

PRECAUTIONS:

The Easy-Flush syringe is designed for use with Luer-compatible connections for intravenous applications. Careful handling is required to avoid damage. Inspect for visible defects or cracks before use. The syringe has graduated markings for precise administration of the saline solution. Do not apply excessive force or use surgical instruments (e.g., forceps) on the product, as this may damage it. Store at a controlled temperature (10 °C to 40 °C) with humidity levels between 55% and 75%. Do not expose to freezing temperatures. Before use, check compatibility with a 0.9% sodium chloride solution. Verify compatibility according to the medication manufacturer's instructions or flush the vascular access device with an appropriate solution to remove medication residues as per manufacturer and institutional guidelines. Healthcare professionals should consider the patient's health status, treatment needs, age, and weight when deciding to flush with a 0.9% sodium chloride solution, particularly in cases of restricted sodium or fluid intake. When prescribing fluids, saline flushes should be included in the total fluid balance to avoid exceeding recommended intake limits. Follow the instructions for use and institutional guidelines for administering flushes.

WARNING:

Read all instructions carefully before use. Do not use if the packaging is open or damaged. For single use only. Opened products must be discarded. Ensure the product/solution remains sterile: the syringe must be intact, and the cap securely attached. There must be no visible leaks that could compromise sterility. Aseptic techniques must be used.

EXPIRATION DATE:

The product has a shelf life of three years from the date of manufacture. Do not use after the expiry date indicated.

STORAGE:

Store at a controlled temperature (10 °C to 40 °C / 50 °F to 104 °F). Do not freeze.

INSTRUCTIONS FOR USE:

1. Check packaging: Do not use if the plastic packaging is damaged or not intact (Figure 1).
2. Open packaging: Carefully tear along the perforation (Figure 2).
3. Inspect the solution: Check for discoloration, cloudiness, sedimentation, or damage to the syringe (Figure 3).
4. Ensure the cap is securely in place. Then, gently press the plunger to release the plunger seal. Improper activation may make handling difficult or cause the plunger to detach. Never pull the plunger back, as this may contaminate the product (Figure 4).
5. Remove air: Remove the syringe cap under aseptic conditions. Hold the syringe upright to gently expel air (Figure 5).
6. Application: Connect the syringe to the access device per institutional protocol and flush (Figure 6). Inject saline slowly to avoid excessive pressure. If resistance is encountered, do not apply excessive force—do not pull or bend the plunger sideways (Figure 7). After flushing, carefully remove the syringe from the access device.
7. Disposal: Dispose of the syringe properly after use. Discard any unused parts. Never reuse disposable syringes! (Figure 8).

INTENDED PATIENT POPULATION:

This syringe is intended for patients requiring a 0.9% sodium chloride flush, regardless of age or gender, provided there are no contraindications.

INTENDED USER:

The product is intended for use by trained medical professionals, including physicians, general practitioners, paramedics, and experienced nurses.

PACKAGE CONTENTS:

The Easy-Flush syringe is supplied sterile and packaged in plastic or paper wrapping. One box contains 30 flushing syringes. Each carton contains eight (8) boxes.

SAFETY WARNINGS:

Do not reuse! Reuse may lead to infections or other complications. Do not re-sterilize or reprocess. Reprocessing may compromise sterility, biocompatibility, and functional integrity. SPM Medicare accepts no liability for reused or re-sterilized products. Do not use if the product has been exposed to freezing temperatures. Do not use after the expiry date. Do not use if the packaging or contents are damaged. Ensure storage within the specified temperature and humidity conditions. Small parts may pose a choking hazard. Dispose of used syringes according to hospital protocol. Do not use if the solution is cloudy, discolored, contains sediment, or has visible particles.

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